



Human Leptin

250 Tubes

Cat. # HL-81K

**HUMAN LEPTIN RIA KIT
250 TUBES (Cat. # HL-81K)**

I. Intended Use	2
II. Principles Of Procedure	2
III. Reagents Supplied	3
IV. Storage and Stability	4
V. Reagent Precautions	4
VI. Materials Required But Not Provided	5
VII. Specimen Collection And Storage	5
VIII. Assay Procedure	6
IX. Calculations and Transformations	9
X. Interpretation	10
XI. Normal Fasting Range	10
XII. Assay Characteristics	10
XIII. Quality Controls	15
XIV. Replacement Reagents	15
XV. Ordering Information	16
XVI. References	17

HUMAN LEPTIN RIA KIT 250 TUBES (Cat. # HL-81K)

I. INTENDED USE

Leptin is a signaling factor encoded by the obese gene in adipose tissue. Administration of recombinant leptin decreases food intake, increases energy expenditures and promotes weight loss.^{1,2} This Human Leptin Radioimmunoassay³ has been developed to quantitate Human Leptin in plasma, serum and tissue culture media. It is a completely homologous assay since the antibody was raised against highly purified Human Leptin and both the standard and tracer are prepared with Human Leptin. ***This kit is for research purposes only.***

II. PRINCIPLES OF PROCEDURE

In radioimmunoassay, a fixed concentration of labeled tracer antigen is incubated with a constant dilution of antiserum such that the concentration of antigen binding sites on the antibody is limited, for example, only 50% of the total tracer concentration may be bound by antibody. If unlabeled antigen is added to this system, there is competition between labeled tracer and unlabeled antigen for the limited and constant number of binding sites on the antibody. Thus, the amount of tracer bound to antibody will decrease as the concentration of unlabeled antigen increases. This can be measured after separating antibody-bound from free tracer and counting one or the other, or both fractions. A calibration or standard curve is set up with increasing concentrations of standard unlabeled antigen and from this curve the amount of antigen in unknown samples can be calculated. Thus, the four basic necessities for a radioimmunoassay system are: a specific antiserum to the antigen to be measured, the availability of a radioactive labeled form of the antigen, a method whereby antibody-bound tracer can be separated from the unbound tracer, and finally, an instrument to count radioactivity.

The Millipore Human Leptin assay utilizes ¹²⁵I-labeled Human Leptin and a Human Leptin antiserum to determine the level of Leptin in serum, plasma or tissue culture media by the double antibody/PEG technique.

III. REAGENTS SUPPLIED

Each kit is sufficient to run 250 tubes and contains the following reagents.

A. Assay Buffer

0.05M Phosphosaline pH 7.4 containing 0.025M EDTA, 0.08% Sodium Azide, 1% RIA Grade BSA and 0.05% Triton X-100

Quantity: 40 mL/vial

Preparation: Ready to use

B. Human Leptin Antibody

Rabbit anti-Human Leptin Serum in Assay Buffer

Quantity: 26 mL/vial

Preparation: Ready to use

C. ¹²⁵I-Human Leptin

¹²⁵I-Human Leptin Label, HPLC purified (specific activity 135 μ Ci/ μ g)

Lyophilized for stability. Freshly iodinated label contains <3 μ Ci, (<111 kBq) calibrated to the 1st Monday of each month.

Quantity: 27 mL/vial upon hydration

Preparation: Contents Lyophilized. Hydrate with entire contents of Label

Hydrating Buffer. Allow to sit at room temperature for 30 minutes, with occasional gentle mixing.

D. Label Hydrating Buffer

Assay Buffer containing Normal Rabbit IgG as a carrier. Used to hydrate ¹²⁵I-Human Leptin

Quantity: 27 mL/vial

Preparation: Ready to use

E. Human Leptin Standards

Purified Recombinant Human Leptin in Assay Buffer at the following concentration: 100 ng/mL

Quantity: 2 mL/vial

Preparation: Ready to use

F. Quality Controls 1 & 2

Purified Recombinant Human Leptin in Assay Buffer

Quantity: 1 mL/vial

Preparation: Ready to use

G. Precipitating Reagent

Goat anti-Rabbit IgG Serum, 3% PEG and 0.05% Triton X-100 in 0.05M Phosphosaline, 0.025M EDTA, 0.08% Sodium Azide

Quantity: 260 mL/vial

Preparation: Ready to use; chill to 4 °C.

IV. STORAGE AND STABILITY

Upon receipt, unused kit may be stored between 2 and 8°C for short term storage. For prolonged storage (>2 weeks), freeze unused kit at $\leq -20^{\circ}\text{C}$. Lyophilized components upon hydration should be stored at $\leq -20^{\circ}\text{C}$ immediately after use, or discarded. Avoid multiple (>5) freeze/thaw cycles. Refer to date on bottle for expiration when stored at $\leq -20^{\circ}\text{C}$. Do not mix reagents from different kits unless they have the same lot number and are unopened.

V. REAGENT PRECAUTIONS

A. Radioactive Materials

This radioactive material may be received, acquired, possessed and used only by research personnel or clinical laboratories for in vitro research tests not involving internal or external administration of the material, or the radiation there from, to human beings or animals. Its receipt, acquisition, possession, use and transfer are subject to the regulations of the U. S. Nuclear Regulatory Commission (NRC) or of a State with which the Commission has entered into an agreement for the exercise of regulatory authority.

The following are suggested general rules for the safe use of radioactive material. The customer's Radiation Safety Officer is ultimately responsible for the safe handling and use of radioactive material.

1. Wear appropriate personal devices at all times while in areas where radioactive materials are used or stored.
2. Wear laboratory coats, disposable gloves, and other protective clothing at all times.
3. Monitor hands, shoes, and clothing and immediate area surrounding the work station for contamination after each procedure and before leaving the area.
4. Do not eat, drink, or smoke in any area where radioactive materials are stored or used.
5. Never pipette radioactive material by mouth.
6. Dispose of radioactive waste in accordance with NRC rules and regulations.
7. Avoid contaminating objects such as telephones, light switches, doorknobs, etc.
8. Use absorbent pads for containing and easily disposing of small amounts of contamination.
9. Wipe up all spills immediately and thoroughly and dispose of the contaminated materials as radioactive waste. Inform Radiation Safety Officer.

V. REAGENT PRECAUTIONS (continued)

B. Sodium Azide

Sodium Azide has been added to all reagents as a preservative at a concentration of 0.08%. Although it is at a minimum concentration, Sodium Azide may react with lead and copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide build up.

VI. MATERIALS REQUIRED BUT NOT PROVIDED

1. Borosilicate glass tubes, 12 x 75 mm. (NOTE: Polypropylene or polystyrene tubes may be used if the investigator finds that the pellet formation is acceptably stable in their system.)
2. 100 μ L pipet with disposable tips
3. 100 μ L & 1.0 mL repeating dispenser
4. Refrigerated swing bucket centrifuge capable of developing 2,000 - 3,000xg. (Use of fixed-angle buckets is not recommended.)
5. Absorbent paper
6. Vortex mixer
7. Refrigerator
8. Gamma Counter

VII. SPECIMEN COLLECTION AND STORAGE

1. A maximum of 100 μ L per assay tube of serum or plasma can be used, although, 50 μ L per assay tube is adequate for most applications. Tissue culture and other media may also be used.
2. Care must be taken when using heparin as an anticoagulant, since an excess will provide falsely high values.⁴ Use no more than 10 IU heparin per mL of blood collected.
3. Specimens can be stored at 4 °C if they will be tested within 24 hours of collection. For longer storage, specimens should be stored at \leq -20 °C. Avoid multiple (>5) freeze/thaw cycles.
4. Avoid using samples with gross hemolysis or lipemia.

VIII. ASSAY PROCEDURE

Standard Preparation

Use care in opening the Standard vial.

Label seven glass tubes 1, 2, 3, 4, 5, 6 and 7. Add 1.0 mL Assay Buffer to each of the seven tubes. Prepare serial dilutions by adding 1.0 mL of the 100 ng/mL standard to tube 1, mix well and transfer 1.0 mL of tube 1 to tube 2, mix well and transfer 1.0 mL of tube 2 to tube 3, mix well and transfer 1.0 mL of tube 3 to tube 4, mix well and transfer 1.0 mL of tube 4 to tube 5, mix well and transfer 1.0 mL of tube 5 to tube 6, mix well and transfer 1.0 mL of tube 6 to tube 7, mix well.

Note: Do not use a Repeater pipette. Change tip for every dilution. Wet tip with Standard before dispensing. Unused portions of standard should be stored at $\leq -20^{\circ}\text{C}$. Avoid multiple freeze/thaw cycles.

Tube #	Standard Concentration	Volume of Assay Buffer to Add	Volume of Standard to Add
1	50 ng/mL	1.0 mL	1.0 mL of 100 ng/mL
2	25 ng/mL	1.0 mL	1.0 mL of 50 ng/mL
3	12.5 ng/mL	1.0 mL	1.0 mL of 25 ng/mL
4	6.25 ng/mL	1.0 mL	1.0 mL of 12.5 ng/mL
5	3.125 ng/mL	1.0 mL	1.0 mL of 6.25 ng/mL
6	1.56 ng/mL	1.0 mL	1.0 mL of 3.125 ng/mL
7	0.78 ng/mL	1.0 mL	1.0 mL of 1.56 ng/mL

For optimal results, accurate pipetting and adherence to the protocol are recommended.

A. Assay Set-Up, Day One

1. Pipet 300 μL of Assay Buffer to the Non-Specific Binding (NSB) tubes (3-4), 200 μL to Reference (Bo) tubes (5-6), and 100 μL to tubes 7 through the end of the assay.
2. Pipet 100 μL of Standards and Quality Controls in duplicate (see flow chart).
3. Pipet 100 μL of each sample in duplicate.

(NOTE: Smaller volumes of sample may be used when Leptin concentrations are anticipated to be elevated or when sample size is limited. Additional Assay Buffer should be added to compensate for the difference so that the volume is equivalent to 100 μL , e.g., when using 50 μL of sample, add 50 μL of Assay Buffer). Refer to Section IX for calculation modification.
4. Pipet 100 μL of ^{125}I -Human Leptin to all tubes. Important: For preparation, see Section III Part C.
5. Pipet 100 μL of Human Leptin antibody to all tubes except Total Count tubes (1-2) and NSB tubes (3-4).
6. Vortex, cover, and incubate overnight (20-24 hours) at 4°C .

B. Day Two

7. Add 1.0 mL of cold (4 °C) Precipitating Reagent to all tubes (except Total Count tubes).
8. Vortex and incubate 20 minutes at 4 °C.
9. Centrifuge, 4 °C, all tubes [except Total Count tubes (1-2)] for 20 minutes at 2,000-3,000xg. NOTE: If less than 2,000xg is used or if slipped pellets have been observed in previous runs, the time of centrifugation must be increased to obtain a firm pellet (e.g., 40 minutes). Multiple centrifuge runs within an assay must be consistent.

Conversion of rpm to xg:

$$xg = (1.12 \times 10^{-5}) (r) (rpm)^2$$

r = radial distance in cm (from axis of rotation to the bottom of the tube)

rpm = rotational velocity of the rotor

10. Immediately decant the supernate of all tubes except Total Count tubes (1-2), drain tubes for at least 15-60 seconds (be consistent between racks), and blot excess liquid from lip of tubes. NOTE: Invert tubes only one time. Pellets are fragile and slipping may occur.
11. Count all tubes in a gamma counter for 1 minute. Calculate the ng/mL of Human Leptin in unknown samples using automated data reduction procedures (see Section IX).

Assay Flow Chart

Day One						Day Two		
	Step 1	Step 2-3	Step 4	Step 5	Step 6	Step 7	Step 8	Step 9-11
Tube #	Add Assay Buffer	Add Standard / QC/ Sample	Add ¹²⁵ I-Leptin Tracer	Add Leptin Antibody	Vortex, Cover, and Incubate 20-24 hrs at 4°C	Add Precipitating Reagent	Vortex, and Incubate 20 min. at 4°C	Centrifuge for 20 min., Decant, and Count pellets
1,2	----	----	100 µL	----		----		
3,4	300 µL	----	100 µL	----		1.0 mL		
5,6	200 µL	----	100 µL	100 µL		1.0 mL		
7,8	100 µL	100 µL of 0.78 ng/mL	100 µL	100 µL		1.0 mL		
9,10	100 µL	100 µL of 1.56 ng/mL	100 µL	100 µL		1.0 mL		
11,12	100 µL	100 µL of 3.125 ng/mL	100 µL	100 µL		1.0 mL		
13,14	100 µL	100 µL of 6.25 ng/mL	100 µL	100 µL		1.0 mL		
15,16	100 µL	100 µL of 12.5 ng/mL	100 µL	100 µL		1.0 mL		
17,18	100 µL	100 µL of 25 ng/mL	100 µL	100 µL		1.0 mL		
19,20	100 µL	100 µL of 50 ng/mL	100 µL	100 µL		1.0 mL		
21,22	100 µL	100 µL of 100 ng/mL	100 µL	100 µL		1.0 mL		
23,24	100 µL	100 µL of QC 1	100 µL	100 µL		1.0 mL		
25,26	100 µL	100 µL of QC 2	100 µL	100 µL		1.0 mL		
27,28	100 µL	100 µL of unknown	100 µL	100 µL	1.0 mL			
29-n	100 µL	100 µL of unknown	100 µL	100 µL	1.0 mL			

IX. CALCULATIONS AND TRANSFORMATIONS

A. Explanation

The calculations for Human Leptin can be automatically performed by most gamma counters possessing data reduction capabilities or by independent treatment of the raw data using a commercially available software package⁵. Choose weighted 4-parameter or weighted log/logit for the mathematical treatment of the data.

NOTE: Be certain the procedure used subtracts the NSB counts from each average count, except Total Counts, prior to final data reduction.

B. Manual Calculation

1. Average duplicate counts for Total Count tubes (1-2), NSB tubes (3-4), Total Binding tubes (reference, Bo) (5-6), and all duplicate tubes for standards and samples to the end of the assay.
2. Subtract the average NSB counts from each average count (except for Total Counts). These counts are used in the following calculations.
3. Calculate the percentage of tracer bound
(Total Binding Counts/Total Counts) X 100
This should be 35-50%.
4. Calculate the percentage of total binding (%B/Bo) for each standard and sample
 $\%B/Bo = (\text{Sample or Standard/Total Binding}) \times 100$
5. Plot the % B/Bo for each standard on the y-axis and the known concentration of the standard on the x-axis using log-log graph paper.
6. Construct the reference curve by joining the points with a smooth curve.
7. Determine the ng/mL of Human Leptin in the unknown samples and controls by interpolation of the reference curve.

NOTE: When sample volumes assayed differ from 100 μ L, an appropriate mathematical adjustment must be made to accommodate for the dilution factor (e.g., if 50 μ L of sample is used, then calculated data must be multiplied by 2).

X. INTERPRETATION

A. Acceptance Criteria

1. The run will be considered accepted when all Quality Control values fall within the calculated Quality Control Range; if any QC's fall outside the control range, review results with supervisor.
2. If the difference between duplicate results of a sample is >10% CV, repeat the sample.
3. The limit of sensitivity for the Human Leptin assay is 0.437 ng/mL + 2 SD (100 μ L sample size).
4. The limit of linearity for the Human Leptin assay is 100 ng/mL (100 μ L sample size). Any result greater than 100 ng/mL should be repeated on dilution using Assay Buffer as a diluent.

XI. NORMAL FASTING RANGE³

Leptin levels are directly correlated with degree of adiposity.

Mean Leptin Values (BMI ranges 18-25):

Lean Men	3.8 \pm 1.8 μ g/L
Lean Women	7.4 \pm 3.7 μ g/L

Levels rise approximately 2.5 times faster in women per unit BMI as compared to men.³

XII. ASSAY CHARACTERISTICS

A. Sensitivity

The lowest level of Leptin that can be detected by this assay is 0.437 ng/mL + 2 SD when using a 100 μ L sample size.

B. Performance

The following parameters of assay performance are expressed as Mean \pm Standard Deviation.

$$\begin{aligned}ED_{80} &= 1.5 \text{ ng/mL} \pm 0.25 \text{ ng/mL} \\ED_{50} &= 7.1 \text{ ng/mL} \pm 0.76 \text{ ng/mL} \\ED_{20} &= 36.4 \text{ ng/mL} \pm 2.19 \text{ ng/mL}\end{aligned}$$

XII. ASSAY CHARACTERISTICS (continued)

C. Specificity

The specificity (also known as selectivity) of an analytical test is its ability to selectively measure the analyte in the presence of other like components in the sample matrix.

Human Leptin	100%
Rat Leptin	<0.2%
Mouse Leptin	<0.2%
Human Insulin	*
Human Proinsulin	*
Rat Insulin	*
Human C-Peptide	*
Glucagon	*
IGF-1	*

*not detectable

D. Precision

Within and Between Assay Variation

Sample No.	Mean ng/mL	Within % CV	Between % CV
1	4.9	8.3	6.2
2	7.2	4.6	5.0
3	10.4	3.9	4.7
4	15.7	4.7	3.0
5	25.6	3.4	3.6

Within and between assay variation was performed on five human serum samples containing varying concentrations of Human Leptin. Data (mean and % CV) shown are from five duplicate determinations of each serum sample in five separate assays.

XII. ASSAY CHARACTERISTICS (continued)

E. Recovery

Spike & Recovery of Leptin in Human Serum

Sample No.	Leptin Added ng/mL	Observed ng/mL	Expected ng/mL	% Recovery
1	0	4.9	-	-
2	2	7.2	6.9	104
3	5	10.4	9.9	105
4	10	15.7	14.9	105
5	20	25.6	24.9	103

Varying concentrations of Human Leptin were added to five human serum samples and the Leptin content was determined by RIA. Mean of the observed levels from five duplicate determinations in five separate assays are shown. Percent recovery was calculated on the observed vs. expected.

F. Linearity

Effect of Serum Dilution

Sample No.	Volume Sampled	Observed ng/mL	Expected ng/mL	% Of Expected
1	100 μ L	45.7	45.7	100
	75 μ L	45.3		99
	50 μ L	45.6		100
	25 μ L	46.1		101
2	100 μ L	31.2	31.2	100
	75 μ L	31.2		100
	50 μ L	31.3		100
	25 μ L	31.0		99
3	100 μ L	13.8	13.8	100
	75 μ L	13.1		95
	50 μ L	12.5		91
	25 μ L	12.1		88
4	100 μ L	9.1	9.1	100
	75 μ L	8.6		95
	50 μ L	8.7		96
	25 μ L	8.4		92

Aliquots of pooled human serum containing varying concentrations of Leptin were analyzed in the volumes indicated. Dilution factors of 1, 1.33, 2 and 4 representing 100 μ L, 75 μ L, 50 μ L and 25 μ L, respectively, were applied in calculating observed concentrations. Mean Leptin levels and percent of expected for five separate assays are shown.

XII. ASSAY CHARACTERISTICS (continued)

G. Example of Assay Results

This data is presented as an example only and should not be used in lieu of a standard curve prepared with each assay.

Tube #	ID	CPM	Ave CPM	Ave Net CPM	% B/Bo	ng/mL
1	Totals	15857	16060			
2	□	16262				
3	NSB	745	723			
4	□	701				
5	Bo	7614	7629	6906		
6	□	7643				
<u>Standards</u>						
7	0.78 ng/mL	7062	7155	6432	0.931	
8		7248				
9	1.56 ng/mL	6383	6428	5705	0.826	
10		6473				
11	3.125 ng/mL	5459	5487	4764	0.690	
12		5515				
13	6.25 ng/mL	4530	4515	3792	0.549	
14		4500				
15	12.5 ng/mL	3354	3354	2631	0.381	
16		3354				
17	25 ng/mL	2373	2449	1726	0.250	
18		2525				
19	50 ng/mL	2008	1983	1260	0.182	
20		1958				
21	100 ng/mL	1592	1585	862	0.125	
		1578				
<u>Controls/Unknown</u>						
23	QC 1	5053	5079	4356	0.631	4.14
24		5105				
25	QC 2	2978	2910	2187	0.317	17.12
26		2842				
27-n	Unknown					

XIII. QUALITY CONTROLS

Good Laboratory Practice requires that Quality Control specimens be run with each standard curve to check the assay performance. Two levels of controls are provided for this purpose. These and any other control materials should be assayed repeatedly to establish mean values and acceptable ranges. Each individual laboratory is responsible for defining their system for quality control decisions and is also responsible for making this system a written part of their laboratory manual. The ranges for Quality Control 1 and 2 are provided on the card insert or can be located at the Millipore website www.millipore.com/bmia.

Recommended batch analysis decision using two controls (Westgard Rule)⁶:

1. When both controls are within ± 2 SD. Decision: Approve batch and release analyte results.
2. When one control is outside ± 2 SD and the second control is within ± 2 SD. Decision: Hold results, check with supervisor. If no obvious source of error is identified by the below mentioned check of systems, the supervisor may decide to release the results.

Technician check of systems:

1. Check for calculation errors
2. Repeat standards and controls
3. Check reagent solutions
4. Check instrument

XIV. REPLACEMENT REAGENTS

Reagents	Cat #
¹²⁵ I-Human Leptin (<3 μ Ci, <111 kBq)	9081
Label Hydrating Buffer (27 mL)	LHB-81
Human Leptin Standards (2 mL each)	8081-K
Human Leptin Antibody (26 mL)	1081-K
Precipitating Reagent (260 mL)	PR-81
Quality Control 1 & 2 (1 mL each)	6081-K
Assay Buffer (40 mL)	AB-PTR

XV. ORDERING INFORMATION

A. To place an order:

For USA Customers:

Please provide the following information to our customer service department to expedite your telephone, fax or mail order:

1. Your name, telephone and/or fax number
2. Customer account number
3. Shipping and billing address
4. Purchase order number
5. Catalog number and description of product
6. Quantity and product size

NOTE: Appropriate license from NRC (or equivalent) must be on file at Millipore before radioactive orders can be shipped.

TELEPHONE ORDERS:

Toll free US 800-MILLIPORE

FAX ORDERS: (636) 441-8050

MAIL ORDERS: Millipore

6 Research Park Drive

St. Charles, Missouri 63304 U.S.A.

For International Customers:

To best serve our international customers, it is Millipore's policy to sell our products through a network of distributors. To place an order or to obtain additional information about Millipore products, please contact your local distributor.

B. Conditions of Sale

All products are for research or manufacturing use only. They are not intended for use in clinical diagnosis or for administration to humans or animals. All products are intended for *in vitro* use only.

C. Material Safety Data Sheets (MSDS)

Material safety data sheets for Millipore products may be ordered by fax or phone. See Section A above for details on ordering.

XVI. REFERENCES

1. Pelleymounter, M.A., et. al. Effects of the obese gene product on body weight regulation in ob/ob mice. *Science*. 269:540-543, 1995.
2. Maffei, M., et. al. Leptin levels in human and rodent: measurement of plasma leptin and ob RNA in obese and weight-reduced subjects. *Nature Med.* Vol. 1, 11:1155-1611, 1995.
3. Ma, Zhongmin, et al. Radioimmunoassay of leptin in human plasma. *Clinical Chemistry*. 42:942-946, June, 1996.
4. Thorell, J.I. *Scand. J. Clin. Lab. Invest.* 31:187, 1973.
5. Feldman, H. and Rodbard, D. "Mathematical Theory of Radioimmunoassay", in: W.D Odell and Doughaday, W.H. (Ed.), Principles of Competitive Protein-Binding Assays. Philadelphia: J.B. Leppincott Company; pp 158-203, 1971.
6. Westgard, J.O., et. al. A multi-rule Shewhart chart for quality control in clinical chemistry. *Clin. Chem.* 27:493-501, 1981.